


Hyperion™

PFOO

without Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- Symmetrical disc design
- Asymmetrical disc design
- Without hub design
- Screw connection

Sizing

- 4 sizes
- Ø 18 to 34 mm


Benefits

- Functional elasticity
- Enhanced shape memory
- Secure positioning
- Anatomical deployment
- Safe attachment/detachment

- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ PFO Occluder is compatible with the Orbit-A Delivery System



- Ⓐ : Left Atrium (LA) disc Ø
- Ⓑ : Right Atrium (RA) disc Ø
- Ⓒ : Connecting Waist Ø
- Ⓓ : Length of Connecting Waist

Ordering Information

Hyperion™ PFOO without Hub

Catalogue #	RA Disc Ø (mm) (B)	LA Disc Ø (mm) (A)	Connecting Waist Ø (mm) (C)	Length of Connecting Waist (mm) (L)	≥ Orbit-A Delivery System
C-PFO-SW 18	18.0	18.0	3.5	6.0	10F
C-PFO-SW 24	24.0	18.0	4.0	7.0	10F
C-PFO-SW 30	28.0	22.0	4.5	7.0	12F
C-PFO-SW 34	34.0	25.0	5.0	7.0	12F

- (A) : Left Atrium (LA) disc Ø
- (B) : Right Atrium (RA) disc Ø
- (C) : Connecting Waist Ø
- (L) : Length of Connecting Waist

