

# Triton™ Aorta AP

## Mechanical Heart Valve

The Sound of Silence



### Product Information

#### Features

- Open pivot
- High durability / rock solid
- Pure pyrolytic carbon
- Optimal axis of suspension
- Tungsten impregnated leaflets
- Titanium strengthening ring

#### Sizing

- 3 sizes
- 16 to 20 mm

#### Benefits

- Minimal noise
- Low thrombogenicity
- Lifetime
- Improved hemodynamic performance
- Excellent imaging
- Optimal leaflet movement

#### Accessories


- Complete sizer set

The open hinging mechanism reduces drastically the noise of the valve movement, thereby improving quality of patient's and his/her family life.

# Ordering Information

## Triton™ Aorta AP Mechanical Heart Valve

Catalogue #	Tissue Anulus Ø (mm)	Orifice Ø (mm)	Orifice Area (mm)	Overall Geographical H (mm)	Opening Angle (°)
CMHV 16AP	16.2	14.8	1.55	8.7	85
CMHV 18AP	18.2	16.8	2.02	9.7	85
CMHV 20AP	20.2	18.8	2.56	10.9	85



Dedicated sizers, accessories and rotatable rings, allowing easy rotation and intra-operative adjustment.



Trademark: Triton is a trademark of Comed BV, The Netherlands. Triton™ is been CE-certified under CE 0344, © Copyrights of Comed BV, The Netherlands - 2016.  
Comed BV, Abe Lenstra Boulevard 36, 8448 JB Heerenveen, The Netherlands, T. +31 (515) 573 399, F. +31 (515) 760 020, info@comedbv.com, www.comedbv.com