

# Anatomy of the Metis™



## 1:1 TRANSMISSION

One-piece construction by Matchlink™ Welding



## STAINLESS STEEL CORE WIRE

135 cm PTFE coated



## EASY GLIDE

The 25 cm spring coil surface is coated with HydroSkin™ offering optimal sliding effect and enhanced trackability



## NITINOL CORE WIRE

Offering push-flex-track



## ACT™ WELDING

Enhancing 1:1 transmission of torque, push and tactile feedback. Extra radiopaque reference



## STAINLESS STEEL SPRING COIL

Offering supple flexibility and high support



## TO SEE IS TO KNOW

High radiopaque Platinum/Nitinal spring coil (3 cm) + ACT™ Welding point 4 cm from tip.



## SHAPE YOUR OWN

Shapeable and re-shapeable guide wire tips. Standard available in Straight Tip (STR) and J-Tip (J).



## YOUR CHOICE

Offering a selection of three-tip-softness. 0.6 g (Floppy), 1.0 g (Supersoft), 2.0 g (Soft)