

loTM

RADIAL ARTERY COMPRESSION TOURNIQUET



The Double Safety Lock provides a perfect balance of safe hemostasis and patient comfort.



WITH DOUBLE SAFETY LOCK

loTM

RADIAL ARTERY COMPRESSION TOURNIQUET

Step-wise deflation provides a safer and more accurate hemostasis.

- The Luer lock prevents air leakage and makes operation easier and safer.



Safety Lock One unlocked



Safety Lock One Locked

- Our proprietary 2ml controlled air release lock design allows for a precise and step-wise deflation resulting in a safer and more comfortable hemostasis for patients.



Safety Lock Two Open



Safety Lock Two Closed
(2ml air evacuation)

- Velcro strip to host the inflation syringe in place
- Clear visualization through balloon allows for easy control of hemostatic process



ORDERING INFORMATION		PACKAGED TEN PER UNIT BOX	
CATALOGUE NUMBER	BANDAGE LENGTH (cm)	Nominal pressure(ml)	Max. inflation volume(ml)
CRACTL20	24.5	15	18
CRACTL20L	27	15	18
CRACTL20XL	29	15	18