

# Prospero™ OTW

## PTA Balloon Dilatation Catheter



LS1611-010.00

Balloon Diameter Ø (mm)	Pressure atm (kPa)												
	4 (405)	5 (507)	6 (608)	7 (709)	8 (810)	9 (912)	10 (1013)	11 (1114)	12 (1216)	13 (1317)	14 (1418)	15 (1020)	16 (1261)
3.00	2.70	2.76	2.89	3.00	3.06	3.09	3.11	3.15	3.19	3.24	3.28	3.33	3.38
3.50	3.44	3.46	3.48	3.50	3.54	3.58	3.62	3.65	3.68	3.71	3.74	3.77	3.81
4.00	3.74	3.83	3.91	4.00	4.05	4.10	4.15	4.20	4.25	4.29	4.33		
4.50	4.30	4.36	4.42	4.50	4.55	4.60	4.66	4.73	4.78	4.82	4.88		
5.00	4.60	4.73	4.86	5.00	5.08	5.17	5.22	5.28	5.33	5.38	5.44		
5.50	5.29	5.35	5.42	5.50	5.54	5.59	5.64	5.70	5.75	5.82	5.89		
6.00	5.62	5.79	5.91	6.00	6.11	6.20	6.26	6.32	6.38	6.43	6.49		
6.50	6.24	6.31	6.40	6.50	6.57	6.65	6.71	6.78	6.83	6.88	6.95		
7.00	6.55	6.73	6.85	7.00	7.04	7.11	7.17	7.23	7.28	7.34	7.39		
8.00	7.50	7.71	7.93	8.00	8.81	8.26	8.35	8.43	8.51	8.60			
9.00	8.40	8.63	8.83	9.00	9.11	9.20	9.28	9.39	9.47	9.58			
10.00	9.50	9.68	9.79	10.00	10.23	10.40	10.54	10.70					
12.00	11.28	11.60	11.75	12.00	12.07	12.20	12.28						

NP: Nominal Pressure

\* Rated Burst Pressure DO NOT EXCEED. (see back)

## Prospero™ Balloon size related to inflation pressure

