


Callisto™ R-Access II

Radial Introducer Kit



Product Information

Features

- Well balanced push/flex
- Kink resistance
- Smooth transition lining

Range

- 4 F to 7 F
- 7 cm sheath length
- 11 cm sheath length
- Colour coded

Puncture

- 22 Ga needle with IV canula

Guide wire

- Polymer jacket with Hydrosurf™
- 0.021"
- 45 cm length
- Straight flexible tip
- Shapeable tip

Quarter-turn-lock


- Dilator fixation
- Easy withdrawl

Hydrosurf™ coating

- Distal 1/3 part of dilator
- Distal 1/3 part of the sheath

High versatility valve

- No leaking
- No air aspiration
- Easy catheter handling



Ordering Information

Callisto™ R-Access II Radial Introducer Kit

| Catalogue # | Ø (F) | Length (cm) | Needle Gauge | Guide Wire Ø (") | Guide Wire Length (cm) | Guide Wire Tip |
|-------------|-------|-------------|--------------|------------------|------------------------|----------------|
| CRSA040721 | 4 F | 7 cm | 22 Ga | 0.021" | 45 cm | Straight |
| CRSA050721 | 5 F | 7 cm | 22 Ga | 0.021" | 45 cm | Straight |
| CRSA060721 | 6 F | 7 cm | 22 Ga | 0.021" | 45 cm | Straight |
| CRSA070721 | 7 F | 7 cm | 22 Ga | 0.021" | 45 cm | Straight |
| CRSA041121 | 4 F | 11 cm | 22 Ga | 0.021" | 45 cm | Straight |
| CRSA051121 | 5 F | 11 cm | 22 Ga | 0.021" | 45 cm | Straight |
| CRSA061121 | 6 F | 11 cm | 22 Ga | 0.021" | 45 cm | Straight |
| CRSA071121 | 7 F | 11 cm | 22 Ga | 0.021" | 45 cm | Straight |

Callisto™ R-Access 24 cm sheath lengths are available as special order. Please contact Comed BV or your local Comed Alliance Group member for information. Minimum order quantity is 10 units.

