

Elara™

Inflation Device



Product Information

Features

- Symmetrical handling
- Ergonomic design
- High pressure tubing
- Analog/fluorescent display
- Screw design
- Rapid trigger release

Measurements

- Maximum pressure: 30 atm
- Stopcock: 3 way 500 psi
- Barrel volume: 20 ml
- Tubing: length 33 cm
pressure 900 psi

Benefits

- For right and left hand users
- Easy piston grip
- Effective pressure transmission
- Visual pressure control
- Rapid inflation
- Rapid/secure deflation

Ordering Information

Catalogue #	Stopcock	Maximum Pressure	Units per Box
C-P-30S	3 way 500 psi	30 atm	5 units

