

Hyperion™

PDAO - Cone Shape

without Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- Cone shape
- Aligned retention skirt
- Without hub design
- Screw connection

Sizing

- 9 sizes
- Ø 04/06 to 20/22 mm


Benefits

- Functional elasticity
- Enhanced shape memory
- Secure positioning
- Anatomical deployment
- Safe attachment/detachment

- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ PDA Occluder is compatible with the Orbit-V Delivery System



(A) : Retention Skirt Ø

(B) : Pulmonary Artery Ø

(C) : Descending Aorta Ø

(L) : Length of Waist

Ordering Information

Hyperion™ PDAO - Cone Shape without Hub

| Catalogue # | Descending Aorta Ø (mm) ③ | Pulmonary Artery Ø (mm) ② | Retention Skirt Ø (mm) ① | Length of Waist (mm) ④ | ≥ Orbit-V Delivery System |
|------------------|------------------------------|------------------------------|-----------------------------|---------------------------|------------------------------|
| C-PDA-SW-II 0406 | 6.0 | 4.0 | 9.0 | 5.0 | 6F |
| C-PDA-SW-II 0608 | 8.0 | 6.0 | 12.0 | 6.0 | 7F |
| C-PDA-SW-II 0810 | 10.0 | 8.0 | 14.0 | 7.0 | 7F |
| C-PDA-SW-II 1012 | 12.0 | 10.0 | 16.0 | 8.0 | 8F |
| C-PDA-SW-II 1214 | 14.0 | 12.0 | 18.0 | 9.0 | 9F |
| C-PDA-SW-II 1416 | 16.0 | 14.0 | 21.0 | 10.0 | 10F |
| C-PDA-SW-II 1618 | 18.0 | 16.0 | 23.0 | 10.5 | 10F |
| C-PDA-SW-II 1820 | 20.0 | 18.0 | 25.0 | 12.0 | 12F |
| C-PDA-SW-II 2022 | 22.0 | 20.0 | 27.0 | 12.0 | 12F |

- ① : Retention Skirt Ø
- ② : Pulmonary Artery Ø
- ③ : Descending Aorta Ø
- ④ : Length of Waist

