

Hyperion™

VSDO - Multi-Fenestrated Membranous with Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- Asymmetrical disc design
- Large LV disc design
- Narrow waist
- With hub design
- Screw connection

Sizing


- 11 sizes
- Ø 4 to 18 mm

Benefits

- Functional elasticity
- Enhanced shape memory
- Secure positioning
- Anatomical deployment
- Coverage of Multi Fenestrations
- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ VSD Occluder is compatible with the Orbit-V Delivery System



- Ⓐ : Left Ventricular (LV) disc Ø
- Ⓑ : Right Ventricular (RV) disc Ø
- Ⓒ : Connecting Waist Ø
- Ⓓ : Length of Connecting Waist

Ordering Information

Hyperion™ VSDO - Multi-Fenestrated Membranous with Hub

Catalogue #	Connecting Waist Ø (mm) ③	Length of Connecting Waist (mm) ④	LV Disc Ø (mm) ①	RV Disc Ø (mm) ②	≥ Orbit-V Delivery System
C-VSD-PME 04	4.0	3.2	12.0	8.0	6F
C-VSD-PME 05	5.0	4.0	13.0	9.0	7F
C-VSD-PME 06	6.0	4.0	14.0	10.0	7F
C-VSD-PME 07	7.0	4.0	15.0	11.0	8F
C-VSD-PME 08	8.0	4.0	16.0	12.0	8F
C-VSD-PME 09	9.0	4.0	17.0	13.0	8F
C-VSD-PME 10	10.0	4.5	18.0	14.0	9F
C-VSD-PME 12	12.0	4.5	20.0	16.0	9F
C-VSD-PME 14	14.0	4.5	22.0	18.0	10F
C-VSD-PME 16	16.0	5.0	24.0	20.0	10F
C-VSD-PME 18	18.0	5.0	26.0	22.0	12F

- ① : Left Ventricular (LV) disc Ø
- ② : Right Ventricular (RV) disc Ø
- ③ : Connecting Waist Ø
- ④ : Length of Connecting Waist

