

Hyperion™

VSDO - Muscular

with Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- Asymmetrical disc design
- With hub design
- Screw connection

Sizing

- 9 sizes
- 7 mm waist length
- Ø 4 to 18 mm


Benefits

- Functional elasticity
- Enhanced shape memory
- Secure positioning
- Anatomical deployment
- Safe attachment/detachment

- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ VSD Occluder is compatible with the Orbit-V Delivery System



- (A) : Left Ventricular (LV) disc Ø
- (B) : Right Ventricular (RV) disc Ø
- (C) : Connecting Waist Ø
- (L) : Length of Connecting Waist

Ordering Information

Hyperion™ VSDO - Muscular with Hub

Catalogue #	Connecting Waist Ø (mm) ③	Length of Connecting Waist (mm) ④	LV Disc Ø (mm) ①	RV Disc Ø (mm) ②	≥ Orbit-V Delivery System
C-VSD-M 06	6.0	7.0	12.0	10.0	7F
C-VSD-M 07	7.0	7.0	13.0	11.0	8F
C-VSD-M 08	8.0	7.0	14.0	12.0	8F
C-VSD-M 09	9.0	7.0	15.0	13.0	8F
C-VSD-M 10	10.0	7.0	16.0	14.0	9F
C-VSD-M 12	12.0	7.0	18.0	15.0	9F
C-VSD-M 14	14.0	7.0	20.0	17.0	10F
C-VSD-M 16	16.0	7.0	22.0	20.0	10F
C-VSD-M 18	18.0	7.0	24.0	22.0	10F

- ① : Left Ventricular (LV) disc Ø
- ② : Right Ventricular (RV) disc Ø
- ③ : Connecting Waist Ø
- ④ : Length of Connecting Waist

