

Hyperion™

Snare



Product Information

Features

- Nitinol shaft
- Nitinol / Gold Tungsten coil loop wire
- Shape memory

Sizing

- L: 129 cm
- Ø Shaft: 0.8 mm / 0.031"
- Ø Snare circle: 15 and 20 mm
- Angle: 90°


Benefits

Easy retrieval of atraumatic foreign bodies in coronary and peripheral cardiovascular system and the extra-cranial neurovascular anatomy


- Functional elasticity
- Anti kinking
- High radiopacity
- ≥ 0.040" compatible

Delivery

Compatible with catheter systems with a minimal ID of 0.040"



Hyperion™ Snare catching a 0.035" guide wire



Ⓛ : Length of shaft

ⓓ : Ø of snare circle

Ⓣ : Ø of snare shaft

Ⓐ : 90°

Ordering Information

Hyperion™ Snare

Catalogue #	CSN-15	CSN-20
Effective Length (cm)	129	129
Ø of Snare Circle (mm)	15	20
Angle °	90	90
Ø Shaft 9 (mm/inch)	0.8 / 0.031"	0.8 / 0.031"

- Ⓛ : Length of shaft
- ⓓ : Ø of snare circle
- Ⓣ : Ø of snare shaft
- Ⓐ : 90°

