

ORDERING INFORMATION

CATALOGUE NUMBER	TYPE	CATHETER SIZE	CATHETER LENGTH (cm)
CCVCBJ1-14-15	Single Lumen	14G	15
CCVCBJ1-14-20	Single Lumen	14G	20
CCVCBJ1-14-30	Single Lumen	14G	30
CCVCBJ1-16-15	Single Lumen	16G	15
CCVCBJ1-16-20	Single Lumen	16G	20
CCVCBJ1-16-30	Single Lumen	16G	30
CCVCBJ1-18-15	Single Lumen	18G	15
CCVCBJ1-18-20	Single Lumen	18G	20
CCVCBJ1-18-30	Single Lumen	18G	30
CCVCBJ1-20-13	Single Lumen	20G	13
CCVCBJ1-20-20	Single Lumen	20G	20
CCVCBJ2-40-05	Double Lumen	4F	5
CCVCBJ2-40-08	Double Lumen	4F	8
CCVCBJ2-40-13	Double Lumen	4F	13
CCVCBJ2-50-08	Double Lumen	5F	8
CCVCBJ2-50-13	Double Lumen	5F	13
CCVCBJ2-50-20	Double Lumen	5F	20
CCVCBJ2-70-15	Double Lumen	7F	15
CCVCBJ2-70-20	Double Lumen	7F	20
CCVCBJ2-70-30	Double Lumen	7F	30
CCVCBJ3-55-08	Triple Lumen	5.5F	8
CCVCBJ3-55-13	Triple Lumen	5.5F	13
CCVCBJ3-70-15	Triple Lumen	7F	15
CCVCBJ3-70-20	Triple Lumen	7F	20
CCVCBJ3-70-30	Triple Lumen	7F	30

Phobos™

CENTRAL VENOUS CATHETER KIT



SAFE AND CONVENIENT CHANNEL FOR CENTRAL VENOUS ACCESS




Phobos™

CENTRAL VENOUS CATHETER KIT

Phobos™ Central Venous Catheter is suitable for short-term vascular access. Multiple channels allow for convenient pressure monitoring, blood sampling, medication and I.V. injection.

Anti-bacterial and leakage-proof design

Proprietary connection design between catheter and extension tubing effectively prevents leakage, avoids blood or medicine residue and bacteria growth.



Biocompatible latex free polyurethane material

Super flexible TPU material offers excellent biocompatibility, kink resistance and radiopacity.

Guidewire introducer syringe


Blue guidewire introducer syringe with hypotube provides easy guidewire insertion.




Cross Section



Single Lumen



Double Lumen



Triple Lumen

	Catheter Size	Distal Lumen Size	Proximal Lumen Size	Middle Lumen Size
Double Lumen	4F	19Ga	20Ga	--
	5F	18Ga	19Ga	--
	7F	14Ga	18Ga	--
Triple Lumen	5.5F	20Ga	22Ga	22Ga
	7F	16Ga	18Ga	18Ga

Wide specification range

A full range of catheter diameters and lengths for various vessel sizes.

Kit contents

Scalpel

Syringe: 5 ml

Tissue Dilator

Introducer Needle

Injection Site Caps

Catheter Fastener

Blue Introducer Syringe: 5ml

Multi-lumen Indwelling Catheter With Extension Line Clamps

Stainless Steel Guidewire (J Tip) With Guidewire Advancer

