

Y-Connectors

Y-Hemostatic Valve Set



Product Information

Types

- Rotating
- Push-Pull

- Push-Pull with tubing
- 20, 25 and 40 cm tubing
- Pressure 250 psi

Torquer compatible with

- 0.014" guide wire
- 0.018" guide wire

Ordering Information

Catalogue #	Type	Tubing	Units per Box
CYHVSNR	Rotating ¹	-	40
CYHVSNP	Push-Pull ²	-	40
CYHVSNP20	Push-Pull ²	20 cm	40
CYHVSNP25	Push-Pull ²	25 cm	40
CYHVSNP40	Push-Pull ²	40 cm	30

Set contains: Y-Connector + torquer + guide wire loader