

Elara™

Inflation Device



Product Information

Features

- Symmetrical handling
- Ergonomic design
- High pressure tubing
- Analog/fluorescent display
- Screw design
- Rapid trigger release

Measurements

- Maximum pressure: 30 atm
- Stopcock: 3 way 500 psi
- Barrel volume: 20 ml
- Tubing: length 33 cm
pressure 900 psi

Benefits

- For right and left hand users
- Easy piston grip
- Effective pressure transmission
- Visual pressure control
- Rapid inflation
- Rapid/secure deflation

Ordering Information

Catalogue #	Stopcock	Maximum Pressure	Units per Box
C-P-30S	3 way 500 psi	30 atm	5 units



Elara™ +

Inflation Device Kit



Product Information

Elara Inflation Device +
Y-Connector Rotatable
or Y-Connector Push-Pull

Elara Inflation Device: Features

- Symmetrical handling
- Ergonomic design
- High pressure tubing
- Analog/fluorescent display
- Screw design
- Rapid trigger release

Measurements

- Maximum pressure: 30 atm
- Stopcock: 3 way 500 psi
- Barrel volume: 20 ml
- Tubing: length 33 cm
pressure 900 psi

Benefits

- For right and left hand users
- Easy piston grip
- Effective pressure transmission
- Visual pressure control
- Rapid inflation
- Rapid/secure deflation

Y-Connector:

- Rotating - 250 psi
- Incl. torque device + loader
- Push-Pull - 250 psi
- Incl. torque device + loader

Ordering Information

Catalogue #	Y-Connector	Units per Box
C-P-30S-YNR	Rotating ¹	1 unit
C-P-30S-YNP	Push/Pull ²	1 unit

Ordering Information for Special Elara™+:

Elara™+ is also available with Push/Pull Y-Connector with tubing in the length 20, 25, 40 and 50 cm.
Minimum quantity for special ordering: 30 units.



Please contact Comed BV or our Comed Alliance Partner concerning delivery times.

