


Hyperion™

Vascular Plug

with Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- Soft elastic design
- With hub design
- Screw connection

Sizing

- 10 sizes
- Ø 4 to 22 mm

Benefits

- Functional elasticity
- Enhanced shape memory
- Secure positioning
- Anatomical deployment
- Safe attachment/detachment

- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ Vascular Plug is compatible with the Orbit-V Delivery System

Ⓧ : Ø

Ⓛ : Length

Ordering Information

Hyperion™ Vascular Plug with Hub

Catalogue #	Waist Ø (mm) $\text{\textcircled{D}}$	Length of Waist (mm) $\text{\textcircled{L}}$	≥ Orbit-V Delivery System
C-PL 04	4.0	5.0	6F
C-PL 06	6.0	6.5	6F
C-PL 08	8.0	7.5	6F
C-PL 10	10.0	8.5	6F
C-PL 12	12.0	9.5	7F
C-PL 14	14.0	10.5	7F
C-PL 16	16.0	11.0	8F
C-PL 18	18.0	11.0	9F
C-PL 20	20.0	12.0	10F
C-PL 22	22.0	13.0	10F

$\text{\textcircled{D}}$: Ø

$\text{\textcircled{L}}$: Length

