

Match Your Size

The balloon material of Oberon™ NC is a non-compliant Nylon, providing a 'Match Your Size' of < 5% size growth.



Balloon Ø	Atm. 4	Atm. 6	Atm. 8	Atm. 10	Atm. 12	Atm. 14	Atm. 16	Atm. 18	Atm. 20	Atm. 22
2.00	1.83	1.86	1.89	1.93	1.97	2.00	2.04	2.08	2.11	2.15
2.50	2.30	2.34	2.38	2.42	2.46	2.50	2.53	2.57	2.61	2.65
2.75	2.55	2.59	2.63	2.67	2.71	2.75	2.78	2.81	2.85	2.89
3.00	2.80	2.84	2.88	2.92	2.96	3.00	3.04	3.07	3.11	3.15
3.25	3.05	3.10	3.14	3.18	3.21	3.25	3.28	3.30	3.33	3.36
3.50	3.25	3.30	3.35	3.40	3.45	3.50	3.54	3.59	3.64	3.68
3.75	3.49	3.55	3.60	3.65	3.70	3.75	3.80	3.84	3.88	3.93
4.00	3.70	3.76	3.83	3.89	3.95	4.00	4.05	4.09	4.14	4.18
4.50	4.20	4.26	4.32	4.38	4.44	4.50	4.55	4.60	4.65	4.70
5.00	4.69	4.76	4.83	4.89	4.95	5.00	5.05	5.10	5.15	5.20

 Nominal Pressure

 Rate Burst Pressure