

Work Your Size

The balloon material of Oberon™ is a semi-compliant Nylon 12, providing a 'Work Your Size' of $\pm 10\%$ size growth.



Balloon Ø	Atm. 2	Atm. 4	Atm. 6	Atm. 8	Atm. 10	Atm. 12	Atm. 14	Atm. 16	Atm. 18	Atm. 20
1,5	1,38	1,41	1,45	1,5	1,55	1,6	1,65	1,7	1,75	1,8
2	1,86	1,91	1,96	2	2,06	2,12	2,18	2,23	2,29	2,35
2,5	2,35	2,4	2,45	2,5	2,56	2,62	2,68	2,75	2,81	2,86
2,75	2,6	2,67	2,72	2,75	2,82	2,88	2,94	3	3,06	3,12
3	2,84	2,89	2,95	3	3,06	3,12	3,18	3,25	3,3	3,36
3,5	3,39	3,44	3,5	3,56	3,62	3,68	3,74	3,81	3,88	3,95

 Nominal Pressure

 Rate Burst Pressure