

# Callisto™ F-Access

## Femoral Introducer Kit



Quarter-turn-screw dilator locking, provides easy and fast system introduction and easy dilator withdrawal

Tubing with stopcock

45 cm stainless steel J-Tip guide wire. 0.035" / 0.038"

Smooth transition of sheath to tapered dilator for stepless accessibility into vasculature

4F


5F

6F

7F

8F

9F



"J-tip" Guidewire

## Product Information

### Features

- Well balanced push/flex
- Kink resistance
- Smooth transition lining

### Range

- 4 F to 9 F
- 11 cm sheath length
- 23 cm sheath length
- Colour coded

### Guide wire

- Stainless Steel
- 0.035" / 0.038"
- 45 / 70 cm length
- J-Tip

### Quarter-turn-lock

- Dilator fixation
- Easy withdrawal

### High versatility valve

- No leaking
- No air aspiration
- Easy catheter handling

# Ordering Information

## Callisto™ F-Access Femoral Introducer Kit

Catalogue #	Ø (F)	Length (cm)	Needle Gauge	Guide Wire Ø (")	Guide Wire Length (cm)	Guide Wire Tip
CFS041135	4 F	11 cm	18 Ga	0.035"	45 cm	J-Tip
CFS051135	5 F	11 cm	18 Ga	0.035"	45 cm	J-Tip
CFS061135	6 F	11 cm	18 Ga	0.035"	45 cm	J-Tip
CFS071135	7 F	11 cm	18 Ga	0.035"	45 cm	J-Tip
CFS081135	8 F	11 cm	18 Ga	0.035"	45 cm	J-Tip
CFS051138	5 F	11 cm	18 Ga	0.038"	45 cm	J-Tip
CFS061138	6 F	11 cm	18 Ga	0.038"	45 cm	J-Tip
CFS071138	7 F	11 cm	18 Ga	0.038"	45 cm	J-Tip
CFS081138	8 F	11 cm	18 Ga	0.038"	45 cm	J-Tip
CFS091138	9 F	11 cm	18 Ga	0.038"	45 cm	J-Tip
CLS062328	6 F	23 cm	18 Ga	0.038"	70 cm	J-Tip
CLS072338	7 F	23 cm	18 Ga	0.038"	70 cm	J-Tip
CLS082338	8 F	23 cm	18 Ga	0.038"	70 cm	J-Tip

