



Callisto™ R-Access I

Radial Introducer Kit



Product Information

Features

- Well balanced push/flex
- Kink resistance
- Smooth transition lining

Range

- 4 F to 7 F
- 7 cm sheath length
- 11 cm sheath length
- Colour coded

Puncture

- 20 Ga needle

Guide wire

- Stainless Steel
- 0.021"
- 45 cm length
- Straight flexible tip
- Shapeable tip

Quarter-turn-lock

- Dilator fixation
- Easy withdrawal

Hydrosurf™ coating

- Distal 1/3 part of dilator
- Distal 1/3 part of sheath

High versatility valve

- No leaking
- No air aspiration
- Easy catheter handling

Ordering Information

Callisto™ R-Access I Radial Introducer Kit

Catalogue #	Ø (F)	Length (cm)	Needle Gauge	Guide Wire Ø (")	Guide Wire Length (cm)	Guide Wire Tip
CRS040721	4 F	7 cm	20 Ga	0.021"	45 cm	Straight
CRS050721	5 F	7 cm	20 Ga	0.021"	45 cm	Straight
CRS060721	6 F	7 cm	20 Ga	0.021"	45 cm	Straight
CRS041121	4 F	11 cm	20 Ga	0.021"	45 cm	Straight
CRS051121	5 F	11 cm	20 Ga	0.021"	45 cm	Straight
CRS061121	6 F	11 cm	20 Ga <td>0.021"</td> <td>45 cm</td> <td>Straight</td>	0.021"	45 cm	Straight
CRS071121	7 F	11 cm	20 Ga	0.021"	45 cm	Straight

Callisto™ R-Access 24 cm sheath lengths are available as special order. Please contact Comed BV or your local Comed Alliance Group member for information. Minimum order quantity is 10 units.

