


Hyperion™

ASDO

with Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- Oxidised short stainless steel hub
- In-disc and in-waist PET patch
- 2 asymmetrical occluder discs
- Screw connection

Sizing

- 26 sizes
- Even sizes:
Ø 6 to 42 mm / small disc size
- Uneven sizes:
Ø 7 to 19 mm / large disc size

Benefits

- Functional elasticity
- Enhanced shape memory
- Easy positioning
- Anatomical deployment
- Safe attachment/detachment

- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ ASD Occluder is compatible with the Orbit-A Delivery System

- Ⓐ : Left Atrium (LA) disc Ø
- Ⓑ : Right Atrium (RA) disc Ø
- Ⓒ : Connecting Waist Ø
- Ⓓ : Length of Connecting Waist

Ordering Information

Hyperion™ ASDO with Hub

EVEN SIZES / SMALL DISCS

Catalogue #	Connecting Waist Ø (mm) (C)	Length of Connecting Waist (mm) (L)	LA Disc Ø (mm) (A)	RA Disc Ø (mm) (B)	≥ Orbit-A Delivery System
C-ASD 06	6.0	3.0	16.0	14.0	8F
C-ASD 08	8.0	3.0	18.0	16.0	8F
C-ASD 10	10.0	3.0	20.0	18.0	9F
C-ASD 12	12.0	3.0	22.0	20.0	9F
C-ASD 14	14.0	3.0	24.0	22.0	9F
C-ASD 16	16.0	4.0	30.0	26.0	10F
C-ASD 18	18.0	4.0	32.0	28.0	10F
C-ASD 20	20.0	4.0	34.0	30.0	10F
C-ASD 22	22.0	4.0	36.0	32.0	10F
C-ASD 24	24.0	4.0	38.0	34.0	12F
C-ASD 26	26.0	4.0	40.0	36.0	12F
C-ASD 28	28.0	4.0	42.0	38.0	12F
C-ASD 30	30.0	4.0	44.0	40.0	12F
C-ASD 32	32.0	4.0	46.0	42.0	12F
C-ASD 34	34.0	4.0	48.0	44.0	14F
C-ASD 36	36.0	4.0	50.0	46.0	14F
C-ASD 38	38.0	4.5	54.0	50.0	14F
C-ASD 40	40.0	4.5	56.0	52.0	14F
C-ASD 42	42.0	4.5	58.0	54.0	14F

Hyperion™ ASDO with Hub

UNEVEN SIZES / LARGE DISC

Catalogue #	Connecting Waist Ø (mm) (C)	Length of Connecting Waist (mm) (L)	LA Disc Ø (mm) (A)	RA Disc Ø (mm) (B)	≥ Orbit-A Delivery System
C-ASD 07	7.0	3.0	21.0	17.0	8F
C-ASD 09	9.0	3.0	23.0	19.0	8F
C-ASD 11	11.0	3.0	25.0	21.0	9F
C-ASD 13	13.0	3.0	27.0	23.0	9F
C-ASD 15	15.0	4.0	29.0	25.0	9F
C-ASD 17	17.0	4.0	31.0	27.0	10F
C-ASD 19	19.0	4.0	33.0	29.0	10F

- (A) : Left Atrium (LA) disc Ø
- (B) : Right Atrium (RA) disc Ø
- (C) : Connecting Waist Ø
- (L) : Length of Connecting Waist

