


Hyperion™

ASDO

without Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- In-disc and in-waist PET patch
- Without hub design
- 2 asymmetrical occluder discs
- Screw connection

Sizing

- 11 sizes
- Ø 6 to 26 mm

Benefits

- Functional elasticity
- Enhanced shape memory
- Easy positioning
- Anatomical deployment
- Safe attachment/detachment

- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ ASD Occluder is compatible with the Orbit-A Delivery System

- Ⓐ : Left Atrium (LA) disc Ø
- Ⓑ : Right Atrium (RA) disc Ø
- Ⓒ : Connecting Waist Ø
- Ⓓ : Length of Connecting Waist

Ordering Information

Hyperion™ ASDO without Hub

EVEN SIZES / SMALL DISCS

Catalogue #	Connecting Waist Ø (mm) ③	Length of Connecting Waist (mm) ④	LA Disc Ø (mm) ①	RA Disc Ø (mm) ②	≥ Orbit-A Delivery System
C-ASD-SW 06	6.0	3.0	16.0	14.0	8F
C-ASD-SW 08	8.0	3.0	18.0	16.0	8F
C-ASD-SW 10	10.0	3.0	20.0	18.0	9F
C-ASD-SW 12	12.0	3.0	22.0	20.0	9F
C-ASD-SW 14	14.0	3.0	24.0	22.0	9F
C-ASD-SW 16	16.0	4.0	30.0	26.0	10F
C-ASD-SW 18	18.0	4.0	32.0	28.0	10F
C-ASD-SW 20	20.0	4.0	34.0	30.0	10F
C-ASD-SW 22	22.0	4.0	36.0	32.0	10F
C-ASD-SW 24	24.0	4.0	38.0	34.0 <td 12F	
C-ASD-SW 26	26.0	4.0	40.0	36.0	12F

- ① : Left Atrium (LA) disc Ø
- ② : Right Atrium (RA) disc Ø
- ③ : Connecting Waist Ø
- ④ : Length of Connecting Waist

