

# Hyperion™

## PDAO - Cone Shape

with Hub



### Product Information

#### Features

- 72 pre-oxidised Nitinol wires
- Cone shape
- Aligned retention skirt
- With hub design
- Screw connection

#### Sizing


- 9 sizes
- Ø 04/06 to 20/22 mm

#### Benefits

- Functional elasticity
- Enhanced shape memory
- Secure positioning
- Anatomical deployment
- Safe attachment/detachment
  
- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

#### Delivery

The Hyperion™ PDA Occluder is compatible with the Orbit-V Delivery System



(A) : Retention Skirt Ø

(B) : Pulmonary Artery Ø

(C) : Descending Aorta Ø

(L) : Length of Waist

# Ordering Information

## Hyperion™ PDAO - Cone Shape with Hub

Catalogue #	Descending Aorta Ø (mm) ③	Pulmonary Artery Ø (mm) ②	Retention Skirt Ø (mm) ①	Length of Waist (mm) ④	≥ Orbit-V Delivery System
C-PDA-II 0406	6.0	4.0	9.0	5.0	6F
C-PDA-II 0608	8.0	6.0	12.0	6.0	7F
C-PDA-II 0810	10.0	8.0	14.0	7.0	7F
C-PDA-II 1012	12.0	10.0	16.0	8.0	8F
C-PDA-II 1214	14.0	12.0	18.0	9.0	9F
C-PDA-II 1416	16.0	14.0	21.0	10.0	10F
C-PDA-II 1618	18.0	16.0	23.0	10.5	10F
C-PDA-II 1820	20.0	18.0	25.0	12.0	12F
C-PDA-II 2022	22.0	20.0	27.0	12.0	12F

- ① : Retention Skirt Ø
- ② : Pulmonary Artery Ø
- ③ : Descending Aorta Ø
- ④ : Length of Waist

