


Hyperion™

VSDO - Perimembranous Zero Edge

with Hub



Product Information

Features

- 72 pre-oxidised Nitinol wires
- No edge on LV disc
- Various waist length
- With hub design
- Screw connection

Sizing

- 10 sizes
- 4 different waist lengths
- Ø 4 to 16 mm

Benefits

- Functional elasticity
- Enhanced shape memory
- Secure positioning
- No valve interference
- No chordae interference
- Anatomical deployment
- Safe attachment/detachment

- High bio-compatibility
- Fast endothelialisation
- Minimised nickel release

Delivery

The Hyperion™ VSD Occluder is compatible with the Orbit-V Delivery System

- Ⓐ : Left Ventricular (LV) disc Ø
- Ⓑ : Right Ventricular (RV) disc Ø
- Ⓒ : Connecting Waist Ø
- Ⓓ : Length of Connecting Waist

Ordering Information

Hyperion™ VSDO - Perimembranous Zero Edge with Hub

Catalogue #	Connecting Waist Ø (mm) ③	Length of Connecting Waist (mm) ④	LV Disc Ø (mm) ①	RV Disc Ø (mm) ②	≥ Orbit-V Delivery System
C-VSD-PMLLD 04	4.0	3.5	10.0	8.0	7F
C-VSD-PMLLD 05	5.0	3.5	11.0	9.0	7F
C-VSD-PMLLD 06	6.0	4.0	12.0	10.0	7F
C-VSD-PMLLD 07	7.0	4.0	13.0	11.0	8F
C-VSD-PMLLD 08	8.0	4.5	14.0	12.0	8F
C-VSD-PMLLD 09	9.0	5.0	15.0	13.0	9F
C-VSD-PMLLD 10	10.0	5.0	18.0	15.0	9F
C-VSD-PMLLD 12	12.0	5.0	22.0	18.0	10F
C-VSD-PMLLD 14	14.0	5.0	24.0	20.0	10F
C-VSD-PMLLD 16	16.0	5.0	26.0	22.0	12F

- ① : Left Ventricular (LV) disc Ø
- ② : Right Ventricular (RV) disc Ø
- ③ : Connecting Waist Ø
- ④ : Length of Connecting Waist

