

Contrast Medium Injection Manifold Kit

Product Information



- All in one kit:
- 1 x angiography control syringe with rotating Luer connector
 - 1 x 500 psi extension tubing
 - 1 x 250 psi manifold with three side ports
 - 2 x infusion sets (IS-G-V and IS-GV3)

Ordering Information

Catalogue Number	Manifold	Control Syringe	High Pressure Ext. Tubing	Infusion Set
CMK-007	CMF253N	CACS112	CHPT90N	IS-G-V IS-G-V3
CMKL-020	CMF253F	CASC112	CHPT120N	IS-G-V IS-G-V3

* Packaged five (5) per unit box

Note:

CMF253N: 250 psi - 3 side ports - On Type Manifold

CMF253F: 250 psi - 3 side ports - On-Type Manifold