

ORDERING INFORMATION

OTW (WORKING LENGTH OF CATHETER IS 800mm)

BALLOON DIAMETER	BALLOON LENGTH (mm)					
	20mm	40mm	60mm	80mm	100mm	120mm
3.0mm	CPTAOW30020	CPTAOW30040	CPTAOW30060	CPTAOW30080	CPTAOW30100	--
3.5mm	--	--	CPTAOW35060	CPTAOW35080	CPTAOW35100	--
4.0mm	CPTAOW40020	CPTAOW40040	CPTAOW40060	CPTAOW40080	CPTAOW40100	--
5.0mm	CPTAOW50020	CPTAOW50040	CPTAOW50060	CPTAOW50080	CPTAOW50100	CPTAOW50120
6.0mm	CPTAOW60020	CPTAOW60040	CPTAOW60060	CPTAOW60080	CPTAOW60100	--
7.0mm	--	CPTAOW70040	CPTAOW70060	CPTAOW70080	--	--
8.0mm	--	CPTAOW80040	CPTAOW80060	--	--	--
9.0mm	--	CPTAOW90040	--	--	--	--
10.0mm	--	CPTAOW10040	--	--	--	--

OTW (WORKING LENGTH OF CATHETER IS 1350mm)

BALLOON DIAMETER	BALLOON LENGTH (mm)				
	40mm	60mm	80mm	100mm	120mm
3.0mm	--	CPTAOW30060L	CPTAOW30080L	CPTAOW30100L	CPTAOW30120L
3.5mm	--	CPTAOW35060L	CPTAOW35080L	CPTAOW35100L	CPTAOW35120L
4.0mm	CPTAOW40040L	--	--	--	--
5.0mm	CPTAOW50040L	CPTAOW50060L	CPTAOW50080L	CPTAOW50100L	--

Additional balloon length available: 20mm, 60mm, 80mm, 100mm and 120mm for all diameters. Balloons with these lengths can be produced on demand at no additional cost, however a minimum quantity of 10 pieces per order applies.

# Prospero<sup>TM</sup>

PTA BALLOON DILATATION CATHETER




ENHANCES PERFORMANCE IN THE MOST TORTUOUS LESIONS




# Prospero<sup>TM</sup>

PTA BALLOON DILATATION CATHETER

- Balloon lengths available up to 120mm
- Strong pushability combined with excellent flexibility
- Large inflation lumen for fast deflation times
- Enhanced Hydrophilic Coating on both catheter shaft and balloon facilitates crossing the most challenging lesions
- Good refolding behavior allowing a safe retrieval, even after multiple inflations



- Two working lengths are available to meet the needs of physicians



## PRODUCT DESCRIPTION

### OVER THE WIRE (OTW)

Material	Nylon 12
Balloon compliance	Semi-compliant
Marker	2 platinum swaged markers
Shaft outer diameter	5F (1.67mm)
Usable catheter length	80cm or 135cm
Lowest crossing profile	0.060" (1.52mm)
Nominal inflation pressure	7atm
Introducer sheath compatibility	5F (for balloon diameters from 3.0 to 6.0mm) 6F (for balloon diameters from 7.0 to 8.0mm) 7F (for balloon diameters from 9.0 to 12.0mm)
Guidewire compatible	0.035" (0.89mm)

